

MICRODRIVE™ DEDICATED HANDPIECES



D30-100
Drill/bur



D32-100
Sagittal saw



D33-100
Oscillating saw



D34-100
Reciprocating saw



D36-100
Wire driver



Extra long bur guard
Order code 2240

Long bur guard
Order code 2230



Handpiece wash cap
Order code 8110

AIR OPERATED HANDPIECE SPECIFICATIONS

Function	Type code	Order code	Speed	Wire capacity/cannulation	Air supply operating pressure	Accessories
Drill	D30-100	219664	0 - 30000 rpm	-	6 - 7 bar	Stericut R series burs
Sagittal saw	D32-100	219774	0 - 25000 cycles/min	-	6 - 7 bar	Stericut S92 saw blades
Oscillating saw	D33-100	219884	0 - 25000 cycles/min	-	6 - 7 bar	Stericut S92 saw blades
Reciprocating saw	D34-100	219994	0 - 18000 cycles/min	-	6 - 7 bar	Stericut S92 saw blades
Wire driver	D36-100	226004	0 - 2000 rpm	0.7 - 1.6 mm	6 - 7 bar	Stericut K wires and guides

See page 10 for 'Wash Assist' air hose system.



MEDICAL EQUIPMENT SPECIALISTS USA

2555 N.W. 102 Ave. # 104 <> Miami <> FL 33172 <> USA
Tel: 786-401-7719 <> Fax: 305-501-4786 <> www.mes-usa.com





deSouther
MEDICAL

AIR HOSES



The Wash Assist co-axial hose system is recommended where machine washing of air hoses is required.

By releasing the probe end from the hose, which is retained on the hose via the rubber strap, the two hose ends can be joined and sealed prior to machine washing.

	Probe	Model	3 metre hose - Order code	4 metre hose - Order code	5 metre hose - Order code
	MA7 Schrader	EPV/MPZ D30 Series	13300 13600	13400 13700	13500 13800
	AO	EPV/MPZ D30 Series	13310 13610	13410 13710	13510 13810
	Draeger	EPV/MPZ D30 Series	13330 13630	13430 13730	13530 13830
	AGA Inst.	EPV/MPZ D30 Series	13320 13620	13420 13720	13520 13830

Other probe connectors are available on request.

STERICUT®



Ask for Stericut® catalogue sheets for De Soutter Medical's wide range of drills, rotary burs, saw blades and K-Wires.

Stericut® disposable precision cutting accessories save time and help cost containment by eliminating the need to clean, sterilise or inspect them for sharpness prior to use.

All accessories are supplied in boxes of five unless otherwise stated.

