

HUVITZ SLIT LAMP HS-5500 / HS-5000



See the difference by looking through the Huvitz Slit Lamp chosen by opinion leaders in the industry: HS-5000 and HS-5500

ILLUMINATION

The light source is a 12-volt, 30-watt high luminance halogen lamp that provides incredible clarity in both the Tower Illumination and Integrated Illumination type slit lamps.



MICROSCOPE

HS-5000 and HS-5500 deliver crisp images and a wide field of view with the global standard Galilean Converging Binocular optical system.

- 12.5x eye pieces / 6x: 38.5mm, 10x: 22.2mm, 16x: 15.2mm, 25x: 10.5mm, 40x: 6.1mm



DESIGN

The slit lamps are offered in the Tower Illumination type (HS-5000) and the Integrated Illumination type (HS-5500) slit lamp models.

Both types are designed in the industry standard type models used for their proven accuracy and reliability.

MAGNIFICATION CONTROL SYSTEM

The five-position drum-style magnification changer provides a wide range of magnification from 6x to 40x easily accessible by rotating the drum.

The design of this system and the uniquely designed Huvitz optic system allows you to offer a more accurate diagnosis and observation to patients without any image distortion in any magnification level.



YELLOW FILTER

A yellow filter is conveniently located near the ocular for effortless insertion of the fluorescein pattern.

With a control lever, any filters are easily inserted.

(Options include cobalt blue, red free, heat absorption, grey, and yellow.)

INTEGRATED CONTROL

The integrated omni style joystick is simple to control.

A trigger button is conveniently mounted on the joystick for easy image and video capture.

Images and videos can be stored simultaneously if the slit lamp is connected to image devices.



SPECIFICATION

		HS-7000	HS-7500	HS-5000	HS-5500
SLIT ILLUMINATION	Slit length(mm)	0.3~14		0.3~12	0.3~14
	Slit width(mm)	0~14 continuous		0~12 continuous	0~14 continuous
	Slit projection	1.167x		1x	1.167x
	Aperture diaphragms	12.5x			
	Filters	Cobalt blue, Red-free, Grey, Heat absorption and Yellow			
	Slit rotation	0°~180° continuous			
PATIENT'S EYE / PRISM	Angle of incidence	0°, 5°, 10°, 15°, 20°	0°~20° continuous	0°, 5°, 10°, 15°, 20°	0°~20° continuous
	Surface working distance	80mm	66mm	80mm	66mm
MICROSCOPE	Microscope	5 steps			
	Type	Galilean converging binocular			
	Magnification	5 position rotating drum			
	Eyepieces	12.5x, 10x		12.5x	
	Total magnifications	6x, 10x, 16x, 25x, 40x			
	Real fields of view (mm)	38.5, 24, 15, 9, 6		38.5, 22.2, 15.2, 10.5, 6.1	
	Interpupillary adjustment	55mm~80mm			
BASE	Vertical movement	28mm			
	Longitudinal movement	78mm			
	Lateral movement	98mm			
	Fine base	10mm			
	Power supply	Free voltage 100-240V AC			
	Voltage frequency	50 / 60Hz			
	Power consumption	70VA			
	Instrument voltage	12V DC			
	Halogen bulb	12V 30W			
	Fixation point bulb	3.4V 20mA			
	DIGITAL CAMERA	Image sensor	1/2" interline CCD		
Image size		up to 1,388 x 1,036 pixels			
Cell size		4.65 μ m x 4.65 μ m			
Resolution depth		8bit or 12bit Raw RGB, YUV 4:2:2			
Transmit method		IEEE 1394 (6pin)			
Transmit speed		400Mbps			
Frame rate		15fps, 7.5fps, 3.75fps			
Lens mount		C-Mount			
Photographing		External trigger or Software trigger			
Dimension		44mm (W) x 29mm (H) x 63mm (D)			
Power consumption		3W (12V DC, from IEEE 1394 cable)			
HIS-5000	CPU	Pentium IV, 2GHz			
	Memory	512 MB (over 1GB recommended)			
	Video card	ATI Radeon 9200 (128MB) or similar			
	System	Microsoft Windows NT (with service pack 6), Windows 2000 (with service pack 4)			
	Transmit method	Windows XP, Microsoft Internet Explorer 6 (Explorer 7 recommended) Standard IEEE 1394a Firewire port OHCI 1.1 compatible (with auxiliary 12V power connector if connected to a laptop computer)			
	Monitor	LCD or CRT (minimum resolution 1,280x1,024 pixels, over 1,600x1,200 pixels recommended)			

OPTION : PC, Table Designs and details can be changed without prior notice for the purposes of improvement.

Huvitz



MEDICAL EQUIPMENT SPECIALISTS USA

2555 N.W. 102 Ave. # 104 <> Miami <> FL 33172 <> USA
Tel: 786-401-7719 <> Fax: 305-501-4786 <> www.mes-usa.com